



- Off-the-shelf anatomic design based on CT scan data
- Corrects temporal hollowing in patients following a pterional approach to the brain
- Anatomic shape allows for deep placement to the temporalis
- Fully porous to promote tissue ingrowth
- Implant bends and holds desired shape without heating
- Screw, or sutures can be placed in any location desired
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile



CAT#	Description	A	B	C
4469	Left	6mm	44mm	43mm
4490	Right	6mm	44mm	43mm

