


- Shaped like a lower eyelid
- Provides internal support in cases of lower eyelid retraction
- Superior strip has smaller pore structure for a smoother edge
- Thin, tissue-like profile for greater eyelid stability
- Fully porous to promote tissue ingrowth and facilitate implant stabilization
- Screw, or sutures can be placed in any location desired
- Easy to trim with scissors or scalpel
- Individually packaged, MR Safe, Sterile



CAT#	Description	A	B	C
4438	Lower Eyelid Spacer w Superior Strip - Large	43mm	22mm	0.45mm

