



- For small to medium orbital floor defects
- Implant provided with a sterile, non-porous transparent template made of the same Su-Por high-density polyethylene
- Available with or without heat bonded non-porous membrane layer
- Fully porous bone-side layer for tissue ingrowth
- Implant bends and holds desired shape without heating
- Smooth template allows for easy repeated removal during fitting without damage to soft tissues
- Outline and contour of the modified template may be used in modification and placement of implant
- Screw, or sutures can be placed in any location desired
- Easy to trim with scissors or scalpel
- MR Safe, Sterile



CAT#	Description	A	B	C
4478	Small - 0.85mm w Template	24mm	24mm	0.85mm
4479	Medium - 0.85mm w Template	30mm	30mm	0.85mm
4480	Large - 1.5mm w Template	35mm	35mm	0.85mm
4481	Small - 1.5mm w Template	24mm	24mm	1.5mm
4482	Medium - 1.5mm w Template	30mm	30mm	1.5mm
4483	Large - 1.5mm w Template	35mm	35mm	1.5mm
4484	Membrane Small w Template	24mm	24mm	1.0mm
4485	Membrane Medium w Template	30mm	30mm	1.0mm
4486	Membrane Large w Template	35mm	35mm	1.0mm

